




NTNU

Norwegian University of
Science and Technology

TMT4320 Nanomaterials **September 28th, 2016**

- TNN5: Characterization of nanomaterials:
XRD, SAXS, EDX, EELS, Nanoindentation




Last lecture

- Summary of:
 - Bottom-up synthesis: Physical/Vapour phase methods
 - Characterization of nanomaterials: TEM, SEM and XPS
- Characterization of Nanomaterials: STM, AFM, XRD and SAXS


Last lecture

- Characterization of Nanomaterials: STM and AFM
 - Near field microscopy techniques:
 - Scanning tunneling microscopy (STM)
 - Tunnel effect
 - Very thin tip electrically powered → tear the electrons by tunnel effect
 - Measuring electric current → Reconstruct where the electrons are




Last lecture

- Characterization of Nanomaterials: STM and AFM
 - Near field microscopy techniques:
 - Atomic Force Microscopy (AFM)
 - Surface sensing technique → Attractive vs Repulsive → Interatomic forces
 - Detection method: Laser beam → Reflexion when cantilever deflects



Last lecture

- Characterization of Nanomaterials: STM and AFM
 - Near field microscopy techniques:
 - Atomic Force Microscopy (AFM)
 - Three types:
 - Contact mode (repulsive forces)
 - Non contact mode (attractive forces)
 - Tapping mode. Non contact mode with larger amplitudes (both repulsive and attractive forces). Larger resolution
 - AFM vs TEM: Convolution effect
Tip shape influence



Today lecture

- Characterization of Nanomaterials: XRD, SAXS, EDX, EELS, FIM, 3DAP, Nanoindentation
 - X-Ray Diffraction (XRD)
 - Small angle X-Ray scattering (SAXS)
 - Energy Dispersive X-Ray Spectroscopy (EDS/EDX)
 - Electron loss energy spectroscopy (EELS)
 - Field Ion Microscope (FIM)
 - Three-dimensional atom probe (3DAP)
 - Nanoindentation


X-ray diffraction

Technique to study the structure, defects and stresses of solids

Beam of x-ray with wavelength from 0.07 to 0.2 nm is diffracted by crystalline specimen according to Bragg's law:

$$\lambda = 2 d \sin \theta$$

λ = X-Ray wavelength
 θ = diffraction angle
 d = interplanar distance



X-ray diffraction

Identify crystalline phases

Structural characteristics (cell parameters, crystallite sizes, defects, etc)

Non-destructive technique

Sample preparation is easy

Cheap

X-ray diffraction

Detection of homogenous and inhomogenous strains due to their dependence on the Bragg angle


Homogeneous or uniform elastic strain shifts the diffraction peaks, without change in peak profile. Shift change in the peaks means change in the lattice constants.

Inhomogeneous strains vary from crystallite to crystallite or even within a single crystallite. As XRD is an averaged information leads to peak broadening.

Peak broadening can also be due to reduction of crystallite size. This can be determined by peak profile analysis → Rietveld refinements

X-ray diffraction

- Grazing Incidence Small Angle X-Ray Scattering (GISAXS)
 - To a first approximation: average height h , average size d , average separation D



X-ray diffraction: instruments

Different sources : Cu, Mo, Co, synchrotron

Different sample holders: air-sensitive, shape, temperature dependent

Different detectors and filters



XRD lab at DMSE (IMT)

D8 Advance : High Temperature

Da Vinci 1

Da Vinci 2

A Unit

D8 Focus

XRD lab at DMSE (IMT)

D8 Advance : High Temperature

- High/Low temperature/pressure
- MRI TC-Wide Range camera
- Vacuum, inert or reactive gas atmospheres
- 190-400 °C with LT stage
- RT-1200 °C with a choice of two radiant heaters
- RT-1400 °C with Pt-Rh or Ta strip heaters
- 0-20 bar gas pressure



XRD lab at DMSE (IMT)

Da Vinci 1

CuK α radiation

2.5° primary and secondary Soller slits

60-90 position sample changer

LynxEye™ SuperSpeed Detector

Variable divergence slit: "V6" means that the divergence slit (the slit between the X-ray source and the sample) opens automatically such that the illuminated length on the sample always remains 6 mm (when using fixed divergence slit, the illuminated length on the sample changes from long to short - this is just geometry).

XRD lab at DMSE (IMT)

Da Vinci 1

- Well suited for:
 - Fast XRD scans for phase identification
 - Samples that can wait in a queue for hours before being measured (the engineer will load the samples into the instrument every morning and afternoon)
 - Energy-discriminating detector makes it useful also for the fluorescing elements (Fe, Co, Mn)
- Not so good for:
 - Samples that change fast in air (then the D8 Focus may be better to use, due to the booking system)
 - Samples that doesn't fit the standard sample holders
 - Measurements where you want to avoid the $K\alpha_2$ contributions (then the A-unit is recommended)

XRD lab at DMSE (IMT)

Da Vinci 1




XRD lab at DMSE (IMT)

Da Vinci 2

MoK α radiation


Capillary sample stage

Typically 3-155 degrees 2θ ($0.46 \text{ \AA}^{-1} < Q < 17.3 \text{ \AA}^{-1}$, $13.6 \text{ \AA} > d > 0.36 \text{ \AA}$)



XRD lab at DMSE (IMT)

- A Unit
 - Bragg-Brentano geometry
 - Data collection from 0.5 to $140^\circ 2\theta$
 - Quartz primary monochromator
 - Cu $K\alpha_1$ radiation
 - Rotating single sample holder with user selectable rotation speeds
- Well suited for:
 - Collection of data for structure solution and/or Rietveld refinement
- Not so good for:
 - Samples containing fluorescing elements (Fe, Co, Mn), depending on the concentration.




XRD lab at DMSE (IMT)

- D8 Focus
 - Bragg-Brentano geometry
 - $\theta-2\theta$ operating mode
 - Cu K α radiation
 - 9 position sample changer
- Well suited for:
 - Booking system - useful if you wish to measure your samples straight after production etc.
- Not so good for:
 - Samples containing cobalt (due to fluorescence). Then the Da Vinci 1 is recommended since its detector has better energy discrimination.



Lay. Double Perov. vs Single Perov.


LnBaM_2O_6 ; Ln = lanthanide, Y; M = Co, Fe, Mn. vs. $\text{Ln}_{0.5}\text{Ba}_{0.5}\text{MO}_3$



LnBaM_2O_6

LDP: A-site cation order

Tetragonal



$\text{Ln}_{0.5}\text{Ba}_{0.5}\text{MO}_3$

SP: A-site cation disorder

Cubic


Why? → Oxygen vacancies

O^{2-} ion and e^- conduction

Ba presence good for H^+ conduction

2
1

Synthesis



EDS/EDX


Electron Dispersive X-Ray Spectroscopy (EDS/EDX)

Analytical technique used for the elemental analysis or chemical characterization

Attached to SEM and TEM microscopes

Limitation on element accuracy → Heavier elements more accurate than lighter elements → Diffraction technique.


Interaction of electrons with matter



Scanning Electron Microscopy SEM

- ❖ Image formation is because of the secondary and back-scattered Electrons
- ❖ Samples are dehydrated and made conductive.

EDX unit



<http://science.howstuffworks.com/scanning-electron-microscope2.htm>

EDS/EDX

Electron Dispersive X-Ray Spectroscopy (EDS/EDX)

Example:

Journal of
Materials Chemistry A



PAPER

[View Article Online](#)
[View Journal](#) | [View Issue](#)

Cite this: *J. Mater. Chem. A*, 2014, 2, 11114

Insight into surface segregation and chromium deposition on $\text{La}_{0.6}\text{Sr}_{0.4}\text{Co}_{0.2}\text{Fe}_{0.8}\text{O}_{3-\delta}$ cathodes of solid oxide fuel cells

Ling Zhao,^a John Drennan,^b Chun Kong,^c Sudath Amarasinghe^d and San Ping Jiang^{*a}

$\text{La}_{0.6}\text{Sr}_{0.4}\text{Co}_{0.2}\text{Fe}_{0.8}\text{O}_{3-\delta}$ (LSCF) perovskite oxide is one of the most important cathode materials in the development of intermediate temperature solid oxide fuel cells (IT-SOFCs), but vulnerable to chromium deposition and poisoning in the presence of gaseous chromium species from the chromia-forming metallic interconnect. Despite extensive studies on Cr deposition on SOFC cathode materials, there is a lack of direct evidence on the surface chemistry and Cr deposition. Here, the fundamental relationship between the surface segregation and Cr deposition of LSCF cathodes is studied on dense LSCF bar samples using a dual beam high resolution focused ion beam (FIB) and a high resolution scanning electron microscope coupled with EDS. FIB-EDS mapping results clearly indicate the segregation of SrO

EDS/EDX

Journal of
Materials Chemistry A

PAPER



View Article Online
View Journal | View Issue

Cite this: *J. Mater. Chem. A*, 2014, 2, 11114

Insight into surface segregation and chromium deposition on $\text{La}_{0.6}\text{Sr}_{0.4}\text{Co}_{0.2}\text{Fe}_{0.8}\text{O}_{3-\delta}$ cathodes of solid oxide fuel cells

Ling Zhao,^a John Drennan,^b Chun Kong,^c Sudath Amarasinghe^d and San Ping Jiang^{*a}




Fig. 2 SEM micrographs of (a and b) the freshly prepared LSCF and (c and d) LSCF samples after being heat treated at 800 °C in the absence of Cr_2O_3 in air for 96 h; (e) enlarged image of the segregated micron-sized particle as circled in (d); and (f) typical EDS spectrum of LSCF grains in (b). Arrows in (d) indicate segregated isolated submicron-sized particles.

EDS/EDX

Cite this: *J. Mater. Chem. A*, 2014, 2, 11114

Insight into surface segregation and chromium deposition on $\text{La}_{0.6}\text{Sr}_{0.4}\text{Co}_{0.2}\text{Fe}_{0.8}\text{O}_{3-\delta}$ cathodes of solid oxide fuel cells

Ling Zhao,^a John Drennan,^b Chun Kong,^c Sudath Amarasinghe^d and San Ping Jiang^{*a}




Fig. 3 SEM micrograph and EDS mapping of the LSCF surface after being heat-treated in the absence of Cr_2O_3 at 800 °C in dry air for 96 h.

EELS


Electron Energy Loss Spectroscopy (EELS)

Technique very close to EDX → Analyze chemical compositions

It works better for lighter elements than EDX

Attached to a TEM unit → Scattered electrons

Interaction of electrons with matter



EELS

Electron Energy Loss Spectroscopy (EELS) @ NTNU



EELS

Electron Energy Loss Spectroscopy (EELS)

Example



Article

pubs.acs.org/cm

Nanoscale Ordering in Oxygen Deficient Quintuple Perovskite $\text{Sm}_{2-\epsilon}\text{Ba}_{3+\epsilon}\text{Fe}_5\text{O}_{15-\delta}$: Implication for Magnetism and Oxygen Stoichiometry

Nadezhda E. Volkova,[†] Oleg I. Lebedev,[‡] Ludmila Ya. Gavrilova,[†] Stuart Turner,[§] Nicolas Gauquelin,[§] Md. Motin Seikh,^{‡,||} Vincent Caignaert,[‡] Vladimir A. Cherepanov,^{*,†} Bernard Raveau,^{*,‡} and Gustaaf Van Tendeloo[§]

EELS

Nanoscale Ordering in Oxygen Deficient Quintuple Perovskite $\text{Sm}_{2-\epsilon}\text{Ba}_{3+\epsilon}\text{Fe}_5\text{O}_{15-\delta}$: Implication for Magnetism and Oxygen Stoichiometry

Nadezhda E. Volkova,[†] Oleg I. Lebedev,[‡] Ludmila Ya. Gavrilova,[†] Stuart Turner,[§] Nicolas Gauquelin,[§] Md. Motin Seikh,^{*,||} Vincent Caignaert,[‡] Vladimir A. Cherepanov,^{*,‡} Bernard Raveau,^{*,‡} and Gustaaf Van Tendeloo[§]




Figure 4. EELS elemental mapping of $\text{Sm}_{2-\epsilon}\text{Ba}_{3+\epsilon}\text{Fe}_5\text{O}_{15-\delta}$. (a) Overview HAADF-STEM image, (b) Fe- $L_{2,3}$ map, (c) Sm- $M_{4,5}$ map, (d) Ba- $M_{4,5}$ map, (e) O- K map, and (f) color overlay with Sm in yellow and Ba in red together with a structural model.

EELS

Nanoscale Ordering in Oxygen Deficient Quintuple Perovskite $\text{Sm}_{2-\epsilon}\text{Ba}_{3+\epsilon}\text{Fe}_5\text{O}_{15-\delta}$: Implication for Magnetism and Oxygen Stoichiometry

Nadezhda E. Volkova,[†] Oleg I. Lebedev,[‡] Ludmila Ya. Gavrilova,[†] Stuart Turner,[§] Nicolas Gauquelin,[§] Md. Motin Seikh,^{‡,||} Vincent Caignart,[‡] Vladimir A. Cherepanov,^{*,‡} Bernard Raveau,^{*,‡} and Gustaaf Van Tendeloo[§]




Figure 6. EELS fine structure of $\text{Sm}_{2-\epsilon}\text{Ba}_{3+\epsilon}\text{Fe}_5\text{O}_{15-\delta}$; (Left panel) O–K edge fine structure signatures from the regions indicated in the central panel. (Center panel) Structural model with indicated EELS integration areas. (Right panel) Fe-L_{2,3} fine structure signatures from the regions indicated in the central panel with references for 4-, 5-, and 6-fold coordinated Fe³⁺, and a simultaneously acquired and energy-shifted Ba M₅ edge.

NanoIndentation

Relatively new technique to obtain mechanical properties of nanometric regions





Fig. 5.15 Schematic diagram of nanoindentation mechanism.



Zr-Cu-Al metallic glass

NanoIndentation

Load-displacement curve




Fig. 5.16 Load-displacement curve of nanoindentation technique (Source: http://commons.wikimedia.org/wiki/File:Load-displacement_curve_%2B_onderdelen.JPG).

NanoIndentation

Zr-Pt alloy




Fig. 5.17 Load–displacement curves for nanocrystalline Zr–Pt alloy (Source: BS Murty, IIT Madras).

NanoIndentation

Example

Journal of Power Sources 273 (2015) 522–529



ELSEVIER

Contents lists available at ScienceDirect

Journal of Power Sources

journal homepage: www.elsevier.com/locate/jpowsour



In-situ Young's moduli of the constitutive layers in a solid oxide fuel cell



Amit Pandey^{a, b, *}, Amit Shyam^a, Zhien Liu^b, Richard Goettler^b

^a Materials Science and Technology Division, Oak Ridge National Laboratory, Oak Ridge, TN 37831, USA

^b LG Fuel Cell Systems Inc., North Canton, OH 44720, USA

H I G H L I G H T S

- A methodology to determine the modulus of functional layers in SOFCs is reported.
- *In-situ* Young's modulus for various functional layers is reported.
- The proposed testing methodology could be applied for other multilayer systems.

NanoIndentation

Example




Fig. 5. Optical micrograph of typical nano indentation imprints on various layers in a substrate supported solid oxide fuel cell (SOFC).

NanoIndentation

Example




Fig. 4. a: Uniaxial monotonic tensile stress-strain response of bilayer specimen i.e. MMA substrate (S) + porous anode barrier (PAB). b: Uniaxial monotonic tensile stress-strain response of trilayer specimen i.e. MMA substrate (S) + porous anode barrier (PAB) + anode current collector (ACC). c: Uniaxial monotonic tensile stress-strain response of multilayer specimen i.e. MMA substrate (S) + porous anode barrier (PAB) + anode current collector (ACC) + anode (A). d: Uniaxial monotonic tensile stress-strain response of multilayer specimen i.e. MMA substrate (S) + porous anode barrier (PAB) + anode current collector (ACC) + anode (A) + electrolyte (E). e: Uniaxial monotonic tensile stress-strain response of multilayer specimen i.e. MMA substrate (S) + porous anode barrier (PAB) + anode current collector (ACC) + anode (A) + electrolyte (E) + cathode (C). f: Uniaxial monotonic tensile stress-strain response of multilayer specimen i.e. MMA substrate (S) + porous anode barrier (PAB) + anode current collector (ACC) + anode (A) + electrolyte (E) + cathode (C) + cathode current collector (CCC).

Summary

- Characterization of Nanomaterials: TEM, SEM, XPS, STM, AFM, XRD, SAXS, EDX, EELS, Nanoindentation
 - Transmission/Scanning Electron Microscopy (TEM/SEM)
 - Scanning tunneling Microscopy
 - Atomic Force Microscopy
 - X-Ray Diffraction (XRD)
 - Small angle X-Ray scattering (SAXS)
 - Energy Dispersive X-Ray Spectroscopy (EDS/EDX)
 - Electron loss energy spectroscopy (EELS)
 - Nanoindentation
 - Many others: FIM, SPM, IR, RAMAN, HR-TEM, ED, SIMS, LEIS, ND, etc...